



E.S.V.I

Vetolabo

Partenaire de vos recherches

Polyethylene (PE) Tubing

Reference	Designation
VET51150	Tubing PE,10 10ft – (.28x.64mm)
VET51151	Tubing PE,10 100ft (.28x.64mm)
VET51154	Tubing PE20, 10ft (.38x1.09mm)
VET51155	Tubing PE20, 100ft (.38x1.09mm)
VET51158	Tubing PE50, 10ft (.58x.99mm)
VET51159	Tubing PE50, 100ft (.58x.99mm)
VET51162	Tubing PE60, 10ft (.76x1.22mm)
VET51163	Tubing PE60, 100ft (.76x1.22mm)
VET51166	Tubing PE90, 10ft (.86x1.32mm)
VET51167	Tubing PE90, 100ft (.86x1.32mm)
VET51170	Tubing PE100, 10ft (.86x1.52mm)
VET51171	Tubing PE10, 10ft (.86x1.52mm)
VET51174	Tubing PE160, 10ft (1.14x1.57mm)
VET51175	Tubing PE160, 100ft (1.14x1.63mm)
VET51178	Tubing PE190, 10ft (1.19x1.70mm)
VET51179	Tubing PE190, 10ft (1.19x1.70mm)

Multipliez la valeur longueur ft par 30,48cm exemple VET51150 (10x30.48) = 3.04m

ESVI VETOLABO

6 Avenue des Andes – Mini-Parc Bâtiment 8 – ZA de Courtabœuf - 91940 Les Ulis cedex

Tél +33 (0)1 60 11 53 44 – Portable +33(0)7 64 36 98 35

91940 Les Ulis cedex

Franck Paolicelli