

*New
Upgrade*

BT100-3J Laboratorial Peristaltic Pump



RoHS

01

High electromagnetic compatibility, reliability and stability

02

Compact structure and simple shape

03

**Easy installation and fixation
Easy cleaning and relocation**

04

Configurable with various pump heads, able to supply a flow at the rate of 0-380 ml/min



BT100-3J - Répond aux exigences de compatibilité électromagnétique de l'UE et est certifié par TUV.

Tête de pompe	Tube	Max débit (ml/min)
YZ1515x – YZII15	13# 14# 19# 25# 17# 18#	380
YZ2515x YZII25	15# 24#	270
DG (1-2) 6 rollers	ID≥3.17mm	48
DG (1-2) 10 rollers	Epaisseur 0.8-1mm	32
DG15-24	16# 25# 17#	300
DG15-28	13# 14# ID≥3.17mm	75
BZ15-13	14# 16# 25# 17#	300
BZ25-13	24#	300
DMD15-13	2*13# 2*14# 2*19# 2*16# 2*25#	350

Débits : 0.0002ml/min-380ml/min

