

CATALOGUE NUMBER	SIZE	COLOUR	G	INTERMEDIATE PACKAGING	BULK PACKAGING
IN-03X12*	0,3mm x 12mm	yellow	30	100 pcs	50 x 100 pcs
IN-03X13	0,3mm x 13mm	yellow	30	100 pcs	50 x 100 pcs
IN-04X19	0,4mm x 19mm	medium grey	27	100 pcs	50 x 100 pcs
IN-045X12	0,45mm x 12mm	brown	26	100 pcs	50 x 100 pcs
IN-045X16	0,45mm x 16mm	brown	26	100 pcs	50 x 100 pcs
IN-045X22	0,45mm x 22mm	brown	26	100 pcs	50 x 100 pcs
IN-05X16	0,5mm x 16mm	orange	25	100 pcs	50 x 100 pcs
IN-05X20*	0,5mm x 20mm	orange	25	100 pcs	50 x 100 pcs
IN-05X25	0,5mm x 25mm	orange	25	100 pcs	50 x 100 pcs
IN-05X40	0,5mm x 40mm	orange	25	100 pcs	50 x 100 pcs
IN-06X25	0,6mm x 25mm	deep blue	23	100 pcs	50 x 100 pcs
IN-06X30	0,6mm x 30mm	deep blue	23	100 pcs	50 x 100 pcs
IN-06X40	0,6mm x 40mm	deep blue	23	100 pcs	50 x 100 pcs
IN-07X30	0,7mm x 30mm	black	22	100 pcs	50 x 100 pcs
IN-07X40	0,7mm x 40mm	black	22	100 pcs	50 x 100 pcs
IN-07X50*	0,7mm x 50mm	black	22	100 pcs	50 x 100 pcs
IN-08X25	0,8mm x 25mm	deep green	21	100 pcs	50 x 100 pcs
IN-08X40	0,8mm x 40mm	deep green	21	100 pcs	50 x 100 pcs
IN-08X50*	0,8mm x 50mm	deep green	21	100 pcs	50 x 100 pcs
IN-09X25	0,9mm x 25mm	yellow	20	100 pcs	50 x 100 pcs
IN-09X40	0,9mm x 40mm	yellow	20	100 pcs	50 x 100 pcs
IN-11X40	1,1mm x 40mm	cream	19	100 pcs	50 x 100 pcs
IN-12X40	1,2mm x 40mm	pink	18	100 pcs	50 x 100 pcs
IN-12X50	1,2mm x 50mm	pink	18	100 pcs	50 x 100 pcs
IN-16X40	1,6mm x 40mm	white	16	100 pcs	50 x 100 pcs
IN-18X40	1,8mm x 40mm	blue-grey	15	100 pcs	50 x 100 pcs
IN-21X40	2,1mm x 40mm	pale green	14	100 pcs	50 x 100 pcs

INDICATION:

- Device designed for intravenous, intramuscular and hypodermic injections

FEATURES:

- Needle tube is made of smooth and durable stainless steel
- Sheath and hub made of polypropylene
- Thin walls assure adequate stream pace
- Tip sharpened in three surfaces (back cut)
- Bevel type: LB/BL standard long bevel
- Small injection and sliding force
- Surface of needle tube covered with sliding agent (silicone)
- Glued needle elements to prevent their separation (glue - epoxy resin)
- Hubs coloured in accordance with ISO code allowing fast recognition of needle size
- Latex-free, phthalate-free
- Non-toxic, non-pyrogenic
- EO sterilized
- Shelf life: 5 years
- Packaging: 1 pc./paper-foil



SIZE AND REFERENCE